

Cable Wiring Guide

1224B - LaserSource Cable, 2A, 2m

Length: 2 meter


Instrument End

Connector: DB-9 Male

Label: Line 1: 1224C
Line 2: LaserSource, 2A, 2m

Pin	Pair	Function / Description	Gauge
1	Tie to #2	Interlock	
2	Tie to #1	Interlock	
3	Shield	Earth Ground	22 AWG
4, 5 & 6	BLK	LD Cathode, PD Cathode	22 AWG
7	WHT	PD Anode	22 AWG
8 & 9	RED	LD Anode	22 AWG
Shell	Shield	Earth Ground	22 AWG

DB-9 Male



Device/Mount End

Connector: 3 terminal TO-can socket

Pin	Color	Function / Description
1	RED	LD Anode
2	BLK	LD Cathode, PD Cathode
3	WHT	PD Anode

