

## Cable Wiring Guide


### 1224C - LaserSource Cable, 2A, 2m

Length: 2 meter

#### Instrument End

Connector: DB-9 Male

Label: Line 1: 1224C  
Line 2: LaserSource, 2A, 2m



Pin	Pair	Function / Description	Gauge
1	Tie to #2	Interlock	
2	Tie to #1	Interlock	
3	Shield	Earth Ground	22 AWG
4, 5 & 7	BLK	LD Cathode, PD Anode	22 AWG
6	WHT	PD Cathode	22 AWG
8 & 9	RED	LD Anode	22 AWG
Shell	Shield	Earth Ground	22 AWG

#### Device/Mount End

Connector: 300-1056 (3 terminal TO-can socket)

Pin	Color	Function / Description
1	RED	LD Anode
2	BLK	LD Cathode, PD Anode
3	WHT	PD Cathode

