

## Cable Wiring Guide

### 1262B-17W2 - TECSOURCE Cable, DB15M to 17W2, 10A, 2m


Length: 2 meter

#### Instrument End

Connector: DB-15 Male

Label: Line 1: 1262B-17W2  
Line 2: TECSOURCE, 10A, 2m

Pin	Pair	Function / Description	Gauge
1	RED	TE+	20 AWG
2	RED/WHT	TE+	20 AWG
3	BLK	TE-	20 AWG
4	BLK/WHT	TE-	20 AWG
5	Shield	Earth Ground	22 AWG
6	Shield	Earth Ground	22 AWG
7	GRN	Sensor+	24 AWG
8	WHT	Sensor-	24 AWG
9	RED/YEL	TE+	20 AWG
10	BLK/YEL	TE-	20 AWG
11	YEL	Fan+	24 AWG
12	BLU	Fan-	24 AWG
13	Short to pin 12	(ID Pin, Male Side Only)	Short to pin 12
14	ORG	R Sens+	24 AWG
15	BRN	R Sens-	24 AWG
Shell	Shield	Earth Ground	



#### Device/Mount End

Connector: 17W2 Female

Pin	Color	Function / Description	Gauge
A1	RED	TE +	20 AWG x 3
	RED/WHT		
	RED/YEL		
A2	BLK	TE -	20 AWG x 3
	BLK/WHT		
	BLK/YEL		
1	n/c	n/c	
2	n/c	n/c	
3	WHT	Sensor -	24 AWG
4	GRN	Sensor +	24 AWG
5	n/c	n/c	
6	n/c	n/c	
7	YEL	Fan +	24 AWG
8	n/c	n/c	
9	n/c	n/c	
10	BRN	Rem Sens -	24 AWG
11	ORG	Rem Sens +	24 AWG
12	n/c	n/c	
13	Shield	Earth Ground	22 AWG
14	n/c	n/c	
15	BLU	Fan -	24 AWG
Shell	Shield	Earth Ground	

