

## Cable Wiring Guide

### 1262B - TECSOURCE Cable, 10A, 2m


Length: 2 meter

#### Instrument End

Connector: DB-15 Male  
 Label: Line 1: 1262B  
 Line 2: TECSOURCE, 10A, 2m

Pin	Pair	Function / Description	Gauge
1	RED	TE+	20 AWG
2	RED/WHT	TE+	20 AWG
3	BLK	TE-	20 AWG
4	BLK/WHT	TE-	20 AWG
5	Shield	Earth Ground	22 AWG
6	Shield	Earth Ground	22 AWG
7	GRN	Sensor+	24 AWG
8	WHT	Sensor-	24 AWG
9	RED/YEL	TE+	20 AWG
10	BLK/YEL	TE-	20 AWG
11	YEL	Fan+	24 AWG
12	BLU	Fan-	24 AWG
13	Short to pin 12	(ID Pin, Male Side Only)	Short to pin 12
14	ORG	R Sens+	24 AWG
15	BRN	R Sens-	24 AWG
Shell	Shield	Earth Ground	22 AWG

DB-15 Male



#### Device/Mount End

Connector: DB-15 Female

Pin	Pair	Function / Description	Gauge
1	RED	TE+	20 AWG
2	RED/WHT	TE+	20 AWG
3	BLK	TE-	20 AWG
4	BLK/WHT	TE-	20 AWG
5	Shield	Earth Ground	22 AWG
6	Shield	Earth Ground	22 AWG
7	GRN	Sensor+	24 AWG
8	WHT	Sensor-	24 AWG
9	RED/YEL	TE+	20 AWG
10	BLK/YEL	TE-	20 AWG
11	YEL	Fan+	24 AWG
12	BLU	Fan-	24 AWG
13	n/c		
14	ORG	R Sens+	24 AWG
15	BRN	R Sens-	24 AWG
Shell	Shield	Earth Ground	22 AWG

DB-15 Female

