

Cable Wiring Guide

1264-DB15 - TECSOURCE Cable, 17W2M to DB15F, 10A, 2m


Length: 2 meter

Instrument End

Connector: 17W2 Male
 Pins: 30A Male Insert
 Hood: D-Sub, 25, Metal

Label: Line 1: 1264-DB15
 Line 2: TECSOURCE, 10A, 2m

<i>Pin</i>	<i>Pair</i>	<i>Function / Description</i>	<i>Gauge</i>
A1	RED	TE+	12 AWG
A2	BLK	TE-	12 AWG
1		n/c	
2		n/c	
3	WHT	Sensor 1-	24 AWG
4	GRN	Sensor 1+	24 AWG
5		n/c	
6		n/c	
7	YEL	Fan+	24 AWG
8		n/c	
9		n/c	
10	BRN/WHT	4-wire Sensor 1-	24 AWG
11	BRN	4-wire Sensor 1+	24 AWG
12		n/c	
13		n/c	
14		n/c	
15	YEL/WHT	Fan-	24 AWG
Shell	Shield	Earth Ground	22 AWG



Device/Mount End

Connector: DB15 Female
 Hood: D-Sub, 15, Metal

<i>Pin</i>	<i>Pair</i>	<i>Function / Description</i>	<i>Gauge</i>
1, 2, 9	RED	TE+	12 AWG
3, 4, 10	BLK	TE-	12 AWG
5	Shield	Earth Ground	22 AWG
6	Shield	Earth Ground	22 AWG
7	GRN	Sensor+	24 AWG
8	WHT	Sensor-	24 AWG
11	YEL	Fan+	24 AWG
12	YEL/WHT	Fan-	24 AWG
13		n/c	
14	BRN	R Sens+	24 AWG
15	BRN/WHT	R Sens-	24 AWG
Shell	Shield	Earth Ground	22 AWG

