

Cable Wiring Guide


1264 - TECSOURCE Cable, 30A, 2m

Length: 2 meter

Instrument End

Connector: 17W2 Male


Label: Line 1: 1264
Line 2: TECSOURCE, 30A, 2m



Pin	Pair	Function / Description	Gauge
A1	RED	TE+	12 AWG
A2	BLK	TE-	12 AWG
1	BLU	TEC On LED+	24 AWG
2	RED/WHT	Remote TEC+	24 AWG
3	WHT	Sensor 1-	24 AWG
4	GRN	Sensor 1+	24 AWG
5	ORG/WHT	AD590-	24 AWG
6	VIO/WHT	Sensor 2-	24 AWG
7	YEL	Fan+	24 AWG
8	BLU/WHT	TEC On LED-	24 AWG
9	BLK/WHT	Remote TEC-	24 AWG
10	BRN/WHT	4-wire Sensor 1-	24 AWG
11	BRN	4-wire Sensor 1+	24 AWG
12	ORG	AD590+	24 AWG
13		n/c	24 AWG
14	VIO	Sensor 2+	24 AWG
15	YEL/WHT	Fan-	24 AWG
Shell	Shield	Earth Ground	22 AWG

Device/Mount End

Connector: 17W2 Female



Pin	Pair	Function / Description	Gauge
A1	RED	TE+	12 AWG
A2	BLK	TE-	12 AWG
1	BLU	TEC On LED+	24 AWG
2	RED/WHT	Remote TEC+	24 AWG
3	WHT	Sensor 1-	24 AWG
4	GRN	Sensor 1+	24 AWG
5	ORG/WHT	AD590-	24 AWG
6	VIO/WHT	Sensor 2-	24 AWG
7	YEL	Fan+	24 AWG
8	BLU/WHT	TEC On LED-	24 AWG
9	BLK/WHT	Remote TEC-	24 AWG
10	BRN/WHT	4-wire Sensor 1-	24 AWG
11	BRN	4-wire Sensor 1+	24 AWG
12	ORG	AD590+	24 AWG
13		n/c	24 AWG
14	VIO	Sensor 2+	24 AWG
15	YEL/WHT	Fan-	24 AWG
Shell	Shield	Earth Ground	22 AWG