

Cable Wiring Guide


1261B - TECSOURCE Cable, 5A, 2m, Pigtailed

Length: 2 meter

Instrument End

Connector: DB-15 Male
Label: Line 1: 1261B
Line 2: TECSOURCE, 5A, 2m

Pin	Pair	Function / Description	Gauge
1	RED	TE+	20 AWG
2	RED/WHT	TE+	20 AWG
3	BLK	TE-	20 AWG
4	BLK/WHT	TE-	20 AWG
5	Shield	Earth Ground	22 AWG
6	Shield	Earth Ground	22 AWG
7	GRN	Sensor+	24 AWG
8	WHT	Sensor-	24 AWG
9	n/c		
10	n/c		
11	YEL	Fan+	24 AWG
12	BLU	Fan-	24 AWG
13	n/c		
14	n/c		
15	n/c		
Shell	Shield	Earth Ground	22 AWG



Device/Mount End

Connector: None, bare wire (pigtail is 3")